


2007-03-29 09:29:00

FIG. 1



COPY OF PAPERS
ORIGINALLY FILED



COPY OF PAPERS
ORIGINALLY FILED

FIG. 2A



40032348-073402





FIG.2B






FIG. 3

I. NVPA (NEUROVASCULAR)

- QUADRATURE / PHASED ARRAY MODE
- CONFIGURATION FILE WILL TURN ON ENTIRE ARRAY
- MAY BE USED FOR ALL TYPE OF BRAIN AND NECK IMAGING, EXCLUDING ECHO PLANAR IMAGING (EPI)
- FOV 46CM (MAXIMUM)

RECOMMENDED USES:

- BRAIN AND OR CERVICAL SPINE LOCALIZERS
- IMAGING OF THE CERVICAL SPINE
- IMAGING OF THE CAROTID ARTERIES
- IMAGING OF THE AORTIC ARCH


COPY OF PAPERS
ORIGINALLY FILED



FIG.4

2. FASTBRN (FAST BRAIN)

- QUADRATURE MODE
 - CONFIGURATION FILE WILL TURN ON ONLY THE HEAD COIL SECTION OF THE ARRAY
 - MAY BE USED FOR ECHO PLANAR IMAGING (EPI)
 - MAY BE USED FOR VASCULAR STUDIES OF THE BRAIN TO REDUCE RECON TIME
 - FOV 24 CM
- RECOMMENDED USES:
- EPI SEQUENCES
 - IMAGING OF THE BRAIN
 - IMAGING OF CIRCLE OF WILLIS



COPY OF PAPERS
ORIGINALLY FILED




FIG. 5

3. HRBRN (HIGH RESOLUTION BRAIN ARRAY)

- QUADRATURE/PHASED ARRAY MODE
- CONFIGURATION FILE WILL TURN ON ONLY THE HEAD COIL SECTION OF THE ARRAY
- USE CONFIGURATION WHEN IMAGING THE BRAIN FOR HIGH RESOLUTION
- COMPATIBLE WITH ALL IMAGING SEQUENCES, EXCEPT EPI
- FOV 24 CM

RECOMMENDED USES:

- HIGH RESOLUTION BRAIN IMAGING
- HIGH RESOLUTION CIRCLE OF WILLIS IMAGING



COPY OF PAPERS
ORIGINALLY FILED

